

E-beam Strategies

- Motivation
- Stage movement strategies
- EBL writing strategies

Motivation

Applications of EBL

- mask fabrication (e.g. chromium on glass)
- direct write (rapid prototyping)
- nano devices in R&D
-


→ different writing strategies required

Recommended Literature:

SPIE HANDBOOK OF MICROLITHOGRAPHY,
MICROMACHINING AND MICROFABRICATION Volume 1:
Microlithography, Chapter 2.1

Stage Movement Strategies


stationary stage



fields

versus

moving stage




stripes

"write-on-the-fly"


EBL Writing Strategies

round (gaussian) beam




versus

shaped beam



vector scan



versus

raster scan




EBL Writing Strategies

EBL Methods

strategy	beam	scan mode	stage
1 (Raith)	gaussian	vector	fixed
2 (Etec)	gaussian	raster	moving
3 (Leica)	shaped		fixed

1st Strategy (Raith)


gaussian beam, vector scan, fixed stage




1st Strategy (**Raith**)

gaussian beam, vector scan, fixed stage

meander mode



line mode



1st Strategy (Raith)


gaussian beam, vector scan, fixed stage

- + fast writing of sparse patterns
(unwritten areas are skipped)
- + easy dose variation from shape to shape
- settling time and hysteresis
have to be calibrated
- overhead time caused by
increased stage settling time

→ Applications: nano lithography, R&D, ...

2nd Strategy (Etec)

gaussian beam, raster scan, moving stage



(e.g. used by MEBES (Etec Systems Inc.))

2nd Strategy (Etec)


gaussian beam, raster scan, moving stage

- + very simple
- + very repeatable → calibration possible
- sparse patterns take as long as dense patterns
- difficult to adjust dose during writing


Applications: mask making, R&D, ...

3rd Strategy (Leica)

shaped beam, moving stage



two fields



3rd Strategy (Leica)

shaped beam, moving stage

+ » 10 x faster than equivalent gaussian beam machines

- extremely complex electron optical column
- complicated calibration routines
- resolution and focus varies with shape size

Applications: mask making,
advanced chip development


Extension: Cell projection
(one of the square apertures is replaced by a more complex shape such as a DRAM cell.)

Summary: 3 Strategies

EBL Methods

strategy	beam	scan mode	stage
1 (Raith)	gaussian	vector	fixed
2 (Etec)	gaussian	raster	moving
3 (Leica)	shaped		fixed

Summary: Strategy used by Raith



Applications:

Nano device fabrication, R&D