

SEM Imaging

- Magnification and Deflection System
- Focussing (contamination dots)
- Imaging


Motivation

To be good in e-beam lithography


means to be good in SEM imaging.

Tip Geometry & Crossover

a) LaB_6 gun



b) thermal field emission



(Gersley, J. Appl. Phys. 65 (3), 914 (1989))

Filament Types


| source type | brightness [A/cm ² /sr] | source size | energy spread [eV] | vacuum [torr] |
|-----------------------------|---------------------------------------|----------------|-----------------------|------------------|
| tungsten thermionic | $\sim 10^5$ | 25 μ m | 2-3 | 10^{-6} |
| LaB ₆ thermionic | $\sim 10^6$ | 10 μ m | 2-3 | 10^{-8} |
| thermal FE (Schottky) | $\sim 10^8$ | 20nm | 0.9 | 10^{-9} |
| cold FE | $\sim 10^9$ | 5nm | 0.22 | 10^{-10} |

SPIE HANDBOOK OF MICROLITHOGRAPHY, MICROMACHINING AND
MICROFABRICATION Volume 1: Microlithography, Chapter 2.2

Beam deflection (E-static/magnetic)


Either magnetic or electrostatic fields can be used to focus electrons just as glass lenses are used to focus rays of light.

Electro-static



Fast deflection

Electro-magnetic




EM lenses
more simple


Beam deflection (E-static/magnetic)

Electron lenses have very poor performance (spherical and chromatic aberrations) compared to light lenses; thus electrons must be kept very close to the axis - small α

Spherical aberration



Chromatic aberration



Beam Deflection E-static or B?

Magnetic deflection

- lower aberrations (as with electrostatic lenses (= condensor in gun))
- but max. frequency is limited by the coil resonances to ~10 MHz.
- max. angle limited by off-axis aberrations to about 1 degree
- Trade off between resolution and field size


on-axis aberrations increase with increasing working distance

\hat{U}





deflection angle and thus off-axis aberrations decrease with increasing working distance

⚠ Compromise field size is much smaller than most chips!

Beam Blanker




Focus / Stigmation Correction




Focussing - Contamination Dots

top view




side view



(Images taken during the acceptance test of Raith 200 at the University of Neuchâtel (Switzerland) 4/98)

Focussing = WD **and** Stigmation



(Image taken at KTH Stockholm (Sweden), see <http://www.nanophys.kth.se/nanophys/facilities/nfl/sem-ebeam.htm>)

Influence of Acceleration Voltage

Low (short penetration depth)

- + Clear surface structures
- + Less damage
- + Less charge up
- + Less edge effect
- Lower resolution

High (large penetration depth)

- + Higher resolution
- Unclear surface structures
- More edge effects
- More sample damage
(heating)

(A guide to Scanning Microscope Observation, Jeol web page 1999)

Influence of Beam Current-Aperture

Low (small aperture)

- + Higher resolution
- + Less damage (heating)
- + Larger depth of focus
- Grainy image

High (large aperture)

- + Smooth image
- + Good signal-to-noise
- Deteriorated resolution
- More damage (heating)
- Lower resolution
- Smaller depth of focus

(A guide to Scanning Microscope Observation, Jeol web page 1999)

Influence of Working Distance

Small


- + Higher resolution
- Smaller depth of focus

Large

- + Larger depth of focus
- Lower resolution

([A guide to Scanning Microscope Observation](#), Jeol web page 1999)

Important rule for SEM imaging



Take always the **LOW** magnification images first!!!